

BIO-REIGN™

BIOACTIVE MATRIX

Bio-Reign™ Bioactive
MOLDABLE Bone Graft Matrix



Our Bioactive Solutions

The evolution of our mineral and collagen composite bone grafts has advanced with the launch of our unique bioactive glass, mineral, and collagen composite bone graft solutions. We have a variety of mineral and collagen composite bone grafts with a wide range of adjustable characteristics. We have expanded even further to offer bioactive moldable bone graft solutions.



Our Composition


Our bioactive composite bone graft matrices are a combination of three components: carbonate apatite anorganic bovine bone mineral, 45S5 bioactive glass, and Type I Collagen. When combined, they provide an optimal scaffold to support the body's natural ability to regenerate new bone.

A Perfect Trio



50%

Carbonate Apatite
Anorganic Bone Mineral



20%

Type I
Collagen



30%

45S5 Bioactive Glass

Bio-Reign™ Bioactive MOLDABLE Bone Graft Matrix


Bio-Reign™ Bioactive Moldable Bone Graft Matrix is composed of carbonate apatite anorganic bone mineral, bioactive glass, and Type I collagen that can be molded to fit the bone defect. It is an osteoconductive, bioactive, porous implant that allows for bony ingrowth across the graft site. The bone graft matrix is slowly resorbed and replaced by new bone tissue during the natural healing process.

Why Bio-Reign™ Bioactive Moldable?

- A Perfect Trio of Components—50% Carbonate Apatite anorganic bone mineral, 30% 45S5 Bioactive Glass, 20% Type I Collagen
- Uniform distribution of bioactive glass and mineral particles throughout the matrix, achieved through our proprietary manufacturing process¹

Why Bio-Reign™ Bioactive Glass Component

- 30% is Optimal: Less is more. Bioactive glass is incorporated into Bio-Reign™ within a suggested critical range of 5-40% for optimal osteoblast growth and calcium phosphate formation in a composite²
- Ideal Particle Range: A narrow particle size distribution limited to 100-300µm to provide a more controlled rate of ion dissolution & surface reactivity, and a more consistent rate of bone bonding & proliferation^{3,4}
- Exemplary Particle Size (100-300µm): Larger sized particles may not fully resorb. Smaller particles may resorb away quickly and impede the upregulation of osteoprogenitor cells^{4,5}




Uniform Particle Distribution

An SEM/EDX Analysis of Bio-Reign™ Bioactive Moldable polished cross sections showing mineral and bioactive glass

- 45S5 BIOACTIVE GLASS PARTICLES
- CARBONATE APATITE ANORGANIC BONE MINERAL
- POROUS TYPE I COLLAGEN MATRIX

Moldable Advantage

- 2 for 1 versatility—Upon hydration, the strip conformation can be used in its original shape or optionally molded into alternative shapes to address the unique contours of each defect
- Combined with either autogenous bone marrow or autograft with saline
- Can also be used with autograft as a bone graft extender
- Moldable, flexible, absorbent, resists migration upon irrigation
- A lengthy 40cc size option unlike any other bioactive moldable bone graft



Almost 2x more absorbent than Vitoss® Bioactive Foam¹ –

	Absorbency (ml/g)
Bio-Reign™ Bioactive Moldable	4.59 ± 0.76
Vitoss® Bioactive Foam	2.70 ± 0.35



Why Carbonate Apatite Bone Mineral?

Optimal **Resorption & Remodeling**^{6,7}

Natural Mineral Structure

Similar to Human Bone Mineral

More **Calcium Phosphate Deposition** than β -TCP⁸


Half the crystallinity than HA, **More Soluble**⁹

Independent studies have shown higher **Osteoclastic & Osteoblastic** activity than β -TCP & HA¹⁰

Why 45S5 Bioactive Glass?

Over **30 Years** of Presence in **Tissue Engineering**^{11,12}

- Favorable Environment for bone regeneration and osteoblast attachment¹³
- Ion Exchange & Release—including soluble tetrahedral silica, which may promote rapid bone formation²
- Cell Proliferation & Differentiation—45S5 Bioactive glass has the ability to stimulate the growth & osteogenic differentiation of human primary osteoblasts¹⁴



Composition of **45S5 Bioactive Glass**

%	Component	
45%	Silicon Dioxide	SiO ₂
24.5%	Calcium Oxide	Ca ₂ O
24.5%	Sodium Oxide	Na ₂ O
6%	Phosphorus Pentoxide	P ₂ O ₅

Why Type I Collagen?

Homologous Molecular Structure to Human Collagen¹⁵

- Highly purified for biocompatibility
- 100% resorbable through normal metabolic pathways¹⁶
- Intrinsic hemostatic properties control minor bleeding^{16,17}
- Well-established long clinical history¹⁶
- Binds proteins and cells and retains biological factors¹⁸
- Single most abundant protein in the human body¹⁹




Five Reasons Why Carbonate Apatite is Superior

1

Optimal **Resorption & Remodeling**^{6,7}


- Not fast like beta-tricalcium phosphate (β -TCP)
- Not slow like hydroxyapatite (HA)
- Ideally, the rate of the bone graft resorption is balanced to the rate of bone remodeling
- Carbonate apatite resorption and remodeling are similar to human bone^{6,7}



2

Natural Mineral Structure Similar to Human Bone Mineral

- Pores provide pathways for cell migration and attachment to lay down new bone
- Carbonate apatite is a better osteoconductive material than HA_2O^{20}




3

More **Calcium Phosphate Deposition** than β -TCP⁸

More calcium phosphate is deposited on the carbonate apatite surface as compared to β -TCP⁸

Osteoblasts prefer attaching to calcium phosphate to lay down new bone



4

Half the Crystallinity than HA. **More Soluble**⁹

Carbonate apatite has half the crystallinity than HA, which enables optimal resorption and remodeling because it more easily resorbs⁹



5

Independent Studies have shown Higher **Osteoclastic & Osteoblastic Activity** than β TCP & HA¹⁰


Osteoclasts break down bone

Carbonate apatite shows higher levels of osteoclastic activity than β -TCP & HA¹⁰


Osteoblasts secrete new bone

Osteoblast proteins are most upregulated with carbonate apatite than β -TCP & HA¹⁰

Higher Osteoclastic Activity



Higher Osteoblastic Activity



Ordering Information



Bio-Reign™ Bioactive Moldable Strips

CATALOG No.	LENGTH	WIDTH	THICKNESS	QUANTITY
BR-05	3.2 cm	2 cm	0.8 cm	5 cc, 1 Strip
BR-10	6.25 cm	2 cm	0.8 cm	10 cc, 1 Strip

References

Bio-Reign™ Bioactive Moldable:

1. Data on file at Royal Biologics, Inc.
2. Gerhardt, L., Boccaccini, A.R. (2010). Bioactive Glass-Ceramic Scaffolds for Bone Tissue Engineering. *Materials*, 3, 3867-3910. Retrieved from <https://doi.org/10.3390/ma3073867>
3. Oonishi, H., Kushitani, S., Yasukawa, E., Iwaki, H., Hench, L.L., Wilson, J., Tsuji, E., Sugihara, T. (1997). Particulate Bioglass Compared With Hydroxyapatite as a Bone Graft Substitute. *Clinical Orthopaedics and Related Research*, 334, 316-325, Lippincott-Raven Publishers, Philadelphia, PA.
4. Schepers, E.J.G., Ducheyne, P. (1997). Bioactive glass particles of narrow size range for the treatment of oral bone defects: a 1-24 month experiment with several materials and particle sizes and size ranges. *Journal of Oral Rehabilitation*, 24, 171-181.
5. Lindfors, N. C., Koski, I., Heikkilä, J. T., Mattila, K. and Aho, A. J. (2010), A prospective randomized 140year followup study of bioactive glass and autogenous bone as bone graft substitutes in benign bone tumors. *J. Biomed. Mater. Res.*, 94B, 157-164. doi:10.1002/jbm.b.31636
6. Matsuura, A., Kubo, T., Doi K., Hayashi, K., Morita, K., Yokota, R., Hayashi, H., Hirata, I., Okazaki, M., Akagawa, Y. (2009). Bone formation ability of carbonate apatite-collagen scaffolds with different carbonate contents. *Dental Materials Journal*, 28(2), 234-242.
7. Ellies, L.G., Carter, J.M., Natiella, J.R., Featherstone, J.D.B., Nelson, D.G.A. (1988). Quantitative analysis of early in vivo tissue response to synthetic apatite implants. *J. of Biomed. Mater. Res.*, 22, 137-148.
8. In vitro data on file at Royal Biologics, Inc.
9. Li, S.T., Chen, H.C., Yuen, D., inventors. (2011). United States patent US 8,980,328. Method of Preparing Porous Carbonate Apatite from Natural Bone.
10. Kanayama, K., Sriarj, W., Shimokawa, H., Ohya, K., Doi, Y., Shibutani, T. 2011. Osteoclast and Osteoblast Activities on Carbonate Apatite Plates in Cell Cultures. *J. Biomaterials*, 26, 435-436.
11. Hench, L. L., & Jones, J. R. (2015). Bioactive Glasses: Frontiers and Challenges. *Frontiers in bioengineering and biotechnology*, 3,194. doi:10.3389/fbioe.2015.00194
12. Hench, L.L. (2013). Chronology of Bioactive Glass Development and Clinical Applications. *New Journal of Glass and Ceramics*, 3(2), 67-73. doi: 10.4236/njgc.2013.32011.
13. Hench, L.L., Polak, J.M., Xynos, I.D., Buttery, L.D.K. (2000). Bioactive Materials to Control Cell Cycle. *Materials Research Innovations*, 3, 313-23. doi: 10.1007/s100190000055
14. Xynos, I.D., Hukkanen, M.V., Batten, J.J., Buttery, L.D.K, Hench, L.L., Polak, J.M. (2000). Bioglass 45S5 stimulates osteoblast turnover and enhances bone formation In vitro: Implications and applications for bone tissue engineering. *Calcif Tissue Int.* 67(4), 321-9.
15. Miller, E.J. (1984). Chemistry of the Collagens and Their Distribution. *Extracellular Matrix Biochemistry*, KA Piez, AH Reddi (eds.). 41-82. Elsevier, New York, NY.
16. Li, S.T. (2000). Biomedical Engineering Handbook, In JD Bronzino (Eds.), *Biologic Biomaterials: Tissue Derived Biomaterials (Collagen)* (1st ed.) 2, 42, 1-23, CRC Press, Boca Raton, FL.
17. Jaffe, R., Deykin, D. (1974). Evidence for a Structural Requirement for the Aggregation of Platelets by Collagen. *The Journal of Clinical Investigation*, 53, 875-883.
18. Geiger, M., Li, R.H., Friess, W. (2003). Collagen sponges for bone regeneration with rhBMP-2. *Science Direct / Elsevier*, 55, 1613-1629. <http://doi.org/10.1016/j.addr.2003.08.010>
19. Ott S. (2003, October 21). Collagen and Bone Matrix. Retrieved from <https://depts.washington.edu/bonebio/ASBMRed/matrix.html>
20. Spence, G., Patel, N., Brooks, R., Rushton, N. (in press). Carbonate substituted hydroxyapatite: Resorption by osteoclasts modifies the osteoblastic response. *Wiley InterScience*. Retrieved from <https://doi.org/10.1002/jbm.a.32083>



BIO-REIGNTM
BIOACTIVE MATRIX

Manufactured for Royal Biologics, Inc.
401 Hackensack Ave, Suite 604
Hackensack NJ 07601
RoyalBiologics.com
Phone: 201-488-1549

Vitoss is a registered trademark of Stryker.

Bio-Reign Bioactive Moldable Bone Graft Matrix is FDA cleared.