


BIOLOGICS

Sports Medicine



Verified Osteoinductive | Highest Quality Standards
Optimized Handling | Sterility Assurance Level of 10^{-6}

Fusion
ORTHOPEDICS

VERIFIED OSTEOINDUCTIVE

- Osteoinductivity is defined by the confirmed presence of chondrocytes, osteocytes, bone marrow, cartilage or new bone using the FDA recognized Athymic Rat Model
- In-vivo osteoinductivity testing is performed AFTER gamma sterilization on every donor lot
- Fusion Biologics are not “Dead Bone” - They are osteoactive with viable, influential BMP’s and other growth factors

OPTIMIZED HANDLING


- OsteoForm’s press fit technology allows for compression of graft by greater than 50% (of its original size) and allows OsteoForm to self-anchor, encouraging maximum bone apposition (re-expansion of graft in irregular shaped gaps) while resisting migration
- Geometric compression ensures graft placement, thus, minimizing additional use of fixation

HIGHEST QUALITY STANDARDS

- OsteoForm precision implants offer a brand new standard for highest quality grafting options
- OsteoForm Shape Memory precision implants deliver consistent performance, flexibility, moldability and shape with uniform porosity
- OsteoForm uses a unique and proprietary process that preserves the inter-connected collagen structures of the trabecular bone and the influential “native” BMP’s and growth factors

STERILITY

- OsteoForm is provided sterile with a sterility assurance level (SAL) of 10^{-6}
- Gamma sterilization is done on every donor lot of OsteoForm
- No other shape memory allograft products offer sterility assurance and verified in-vivo osteoinductivity/bone growth



OSTEOFORM HISTOLOGY: Cell Staining From Rat Ossicle Model

- A New Bone
- B Osteocytes
- C OsteoForm
- D Cartilage
- E Chondrocytes
- F Bone Marrow
- G Osteoblasts

ALLOCELL™ BMA/PRP

PRODUCT	DESCRIPTION
AC-BP-0010-S	AlloCell 100 BMA PRP, 1.0cc (Jar)
AC-BP-0025-S	AlloCell 100 BMA PRP, 2.5cc (Jar)
AC-BP-0050-S	AlloCell 100 BMA PRP, 5.0cc (Jar)
AC-BP-0100-S	AlloCell 100 BMA PRP, 10.0cc (Jar)




OSTEOFORM™ SPONGE BLOCK

PRODUCT	DESCRIPTION
OF-SB-1010-S	OsteoForm Sponge Block, 10mm ³
OF-SB-1212-S	OsteoForm Sponge Block, 12mm ³
OF-SB-1414-S	OsteoForm Sponge Block, 14mm ³



OSTEOFORM™ FLEX STRIP

PRODUCT	DESCRIPTION
OF-FS-2010-S	OsteoForm Flex Strip, 20mm x 10mm x 10mm
OF-FS-2015-S	OsteoForm Flex Strip, 20mm x 15mm x 7mm
OF-FS-2510-S	OsteoForm Flex Strip, 25mm x 10mm x 7mm
OF-FS-2619-S	OsteoForm Flex Strip, 26mm x 19mm x 7mm



OSTEOFORM™ FILLER (SYRINGE)

PRODUCT	DESCRIPTION
OF-FS-0025-S	OsteoForm Filler, 2.5cc Syringe
OF-FS-0050-S	OsteoForm Filler, 5.0cc Syringe
OF-FS-0100-S	OsteoForm Filler, 10.0cc Syringe



OSTEOFORM™ FILLER (JAR)

PRODUCT	DESCRIPTION
OF-FJ-0025-S	OsteoForm Filler, 2.5cc Jar
OF-FJ-0050-S	OsteoForm Filler, 5.0cc Jar
OF-FJ-0100-S	OsteoForm Filler, 10.0cc Jar



OSTEOPURE™ 100 PUTTY

PRODUCT	DESCRIPTION
OP-PS-0010-S	OsteoPure 100 Putty, 1.0cc (Syringe)
OP-PS-0025-S	OsteoPure 100 Putty, 2.5cc (Syringe)
OP-PS-0050-S	OsteoPure 100 Putty, 5.0cc (Syringe)
OP-PS-0100-S	OsteoPure 100 Putty, 10.0cc (Syringe)
OP-PJ-0025-S	OsteoPure 100 Putty, 2.5cc (Jar)
OP-PJ-0050-S	OsteoPure 100 Putty, 5.0cc (Jar)
OP-PJ-0100-S	OsteoPure 100 Putty, 10.0cc (Jar)




OSTEOFORM™ TRICORTICAL BLOCK WEDGE

PRODUCT	DESCRIPTION
OF-TB-1618-S	OsteoForm Tricortical Block Wedge Small, 16-18mm
OF-TB-2022-S	OsteoForm Tricortical Block Wedge Medium, 20-22mm
OF-TB-2426-S	OsteoForm Tricortical Block Wedge Large, 24-26mm



PERICARDIUM™ GRAFTS

PRODUCT	DESCRIPTION
BI-PG-0404-S	Pericardium Graft, 4cm x 4cm
BI-PG-0606-S	Pericardium Graft, 6cm x 6cm
BI-PG-0705-S	Pericardium Graft, 7cm x 5cm
BI-PG-1005-S	Pericardium Graft, 10cm x 5cm



DERMAFORM™ DERMAL GRAFTS

PRODUCT	DESCRIPTION
DF-DG-0054-S	DermaForm Dermal Graft, 5.0cm x 4.0cm x 0.5mm
DF-DG-0075-S	DermaForm Dermal Graft, 7.0cm x 5.0cm x 0.5mm
DF-DG-0105-S	DermaForm Dermal Graft, 10.0cm x 5.0cm x 0.5mm
DF-DG-1054-S	DermaForm Dermal Graft, 5.0cm x 4.0cm x 1.0mm
DF-DG-1075-S	DermaForm Dermal Graft, 7.0cm x 5.0cm x 1.0mm
DF-DG-1105-S	DermaForm Dermal Graft, 10.0cm x 5.0cm x 1.0mm
DF-DG-2054-S	DermaForm Dermal Grafts, 5.0cm x 4.0cm x 2.0mm
DF-DG-2075-S	DermaForm Dermal Graft, 7.0cm x 5.0cm x 2.0mm
DF-DG-2105-S	DermaForm Dermal Graft, 10.0cm x 5.0cm x 2.0mm



SPORTS MEDICINE GRAFTS

PRODUCT	DESCRIPTION
BI-BT-1010-S	Bone Tendon Bone - Whole
BI-BT-1020-S	Bone Tendon Bone - Hemi
BI-BT-1030-S	Bone Tendon Bone - Preshaped
BI-AT-2010-S	Achilles Tendon with Bone Block
BI-AT-2020-S	Achilles Tendon - Preshaped
BI-TT-3010-S	Anterior Tibialis Tendon
BI-TT-3020-S	Posterior Tibialis Tendon
BI-ST-4010-S	Semitendinosus Tendon
BI-GT-5010-S	Gracilis Tendon
BI-PT-6010-S	Peroneus Tendon

