



DYNABRIDGE

SUPERELASTIC NITINOL IMPLANT




WHY DYNABRIDGE?

The proof is in the science




WHY NITINOL?

Dynamic compression is ideal




YOCO AKIN BLADE

The YOCO offers an easier, faster, and reproducible wedge cut with only one pass.



The Statistics



PLACEMENT GUIDE

1 Hallux IPJ Arthrodesis
9mm, 11mm, 13mm

2 Akin Osteotomy
9mm, 11mm

3 First MPJ Arthrodesis
15mm, 18mm, 20mm


4 Chevron Osteotomy
11mm, 13mm, 15mm

5 Base Wedge Osteotomy
15mm, 18mm

6 First Metatarsocuneiform Arthrodesis / Lapidus Procedure
15mm, 18mm, 20mm

7 Naviculocuneiform Arthrodesis
15mm, 18mm, 20mm

8 Talonavicular Arthrodesis
15mm, 18mm, 20mm, 25mm



9 Weil Osteotomy
9mm, 11mm, 13mm

10 Bunionette Correction
9mm, 11mm

11 Lisfranc Fusion
15mm, 18mm, 20mm

12 Jones Fracture
9mm, 11mm

13 Calcaneocuboid Arthrodesis
15mm, 18mm, 20mm


14 Evans Calcaneal Osteotomy
15mm, 18mm, 20mm, 25mm

15 Dwyer Osteotomy
15mm, 18mm, 20mm, 25mm

ORDERING INFO

DYNABRIDGE

SUPERELASTIC NITINOL IMPLANT



Before spin

After spin
Mechanical Deformation



Sterile Kit

SIZING

	Size	Tray Implant	Kit Implant	Single Implant	Bridge	Leg A	Leg B
Style A	9x8x8mm	DB-09-0808-A	DB-09-0808-SA	DB-09-0808-SA1	9mm	8mm	8mm
	11x10x10mm	DB-11-1010-A	DB-11-1010-SA	DB-11-1010-SA1	11mm	10mm	10mm
	13x10x10mm	-	DB-13-1010-SA	DB-13-1010-SA1	13mm	10mm	10mm

Style B	15x15x15mm	-	DB-15-1515-SB	DB-15-1515-SA1	15mm	15mm	15mm
	18x18x18mm	-	DB-18-1818-SB	DB-18-1818-SA1	18mm	18mm	18mm
	20x20x20mm	-	DB-20-2020-SB	DB-20-2020-SA1	20mm	20mm	20mm
	25x20x20mm	-	DB-25-2020-SB	DB-25-2020-SA1	25mm	20mm	20mm

	Size	Tray	Kitted	Seperate
YOCO	Akin (3mm)	CD-OP-1003	CD-OP-1003-S	CD-OP-1003-S1

FUSION ORTHOPEDICS USA, LLC.

6859 E Rembrandt Ave #122
Mesa, Arizona 85212 USA
(800)-403-6876

Fusion

ORTHOPEDICS

FDA Cleared.

All Rights Reserved.
DynaBridge is a Trademark
of Fusion Orthopedics USA, LLC.
SLS008 Rev12122023

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Surgeons should rely on their education and professional training, as the Fusion Instru- ments, implants, and techniques ("Products") augment traditional surgeries. Fusion Orthopedics, and its affiliates ("Fusion"), does not dispense medical advice and while Fusion recommends that surgeons be trained in the use of any particular Fusion Product before using it in surgery, surgeons must still rely on their professional training and expertise. The information presented is intended to demonstrate the breadth of Fusion's Products. A surgeon and patient must always refer to the current written surgical technique guide, fusionorthopedics.com and associated websites, package insert, product label, and instructions for use before considering and using any Fusion Product. Products and associated accessories may not be available in all markets because of Product availability, and Products are subject to the regulatory law and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Products in your area. Patients must solely rely on advice from their surgeon regarding the risks and all steps of a particular surgery with the utilization of Fusion's Products. As with any surgery, individual results and recovery time vary. Recovery time may range from months to years. Known risks associated with surgery are, but not limited to, injury, significant scarring, nerve damage, mal-union, long-term pain, implant pain, non-union, misalignment, loss of motion, bone shortening, infections, swelling, stiffness, bone fracture, and subsequent required amputation. These works, and its contents, are protected by U.S. and International copyright laws. Reproduction and/or distribution of these works, and its contents, without the express written consent of Fusion Orthopedics USA, LLC is prohibited. © 2023 Fusion Orthopedics USA, LLC