

# HAMMERTECH

FIXATION SYSTEM



A TRADEMARK OF  
**Fusion**  
ORTHOPEDICS

# HAMMERTECH

FIXATION SYSTEM

## SYSTEM INFO

**Sterile Kit**  
5 Sizes Available for a Total of 15 Different Kits

**Handle**  
Ergonomic Design

**K-Wires**  
Double-Trocar



**Cannulated Driver**  
Drill until hard stop

**Cannulated Drill Bit**  
Precision Fitment

**Implant**  
Cannulated PEEK, Ti Straight, Ti 10° Angled

### Add-Ons:

Sterile Packed and Sold Separately



**Awl**  
Designed for joint preparation



**Removal Tool**  
For easy and efficient dislodgement of PEEK implant



**Cup & Cone Reamers**  
Designed for increased contact surface



**Sizer Wheel**  
To determine optimal implant fitting




**Implant Singles**  
Additional single-packed implants

### Pull Out Force




## SIZING


### PEEK



### Ti Straight



### Ti 10° Angled



	Size	Drill	K-Wire
Black	Petite	2.4mm	1.0mm
Gold	Extra Small	2.4mm	1.0mm
Purple	Small	2.75mm	1.4mm
Green	Medium	2.95mm	1.4mm
Blue	Large	3.25mm	1.6mm

	Size	Drill	K-Wire
Black	Petite	2.3mm	1.0mm
Gold	Extra Small	2.4mm	1.0mm
Purple	Small	2.8mm	1.4mm
Green	Medium	3.0mm	1.4mm
Blue	Large	3.2mm	1.6mm

	Size	Drill	K-Wire
Black	Petite	2.3mm	N/A
Gold	Extra Small	2.4mm	N/A
Purple	Small	2.8mm	N/A
Green	Medium	3.0mm	N/A
Blue	Large	3.2mm	N/A

All dimensions shown are approximate and in mm

**FUSION ORTHOPEDICS USA, LLC.**  
6859 E Rembrandt Ave #122  
Mesa, Arizona 85212 USA  
(800)-403-6876

**Fusion**  
ORTHOPEDICS

**FDA Cleared.**  
All Rights Reserved.  
HammerTech is a Trademark  
of Fusion Orthopedics USA, LLC.  
SLS009 10032023

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Surgeons should rely on their education and professional training, as the Fusion Instruments, implants, and techniques ("Products") augment traditional surgeries. Fusion Orthopedics, and its affiliates ("Fusion"), does not dispense medical advice and while Fusion recommends that surgeons be trained in the use of any particular Fusion Product before using it in surgery, surgeons must still rely on their professional training and expertise. The information presented is intended to demonstrate the breadth of Fusion's Products. A surgeon and patient must always refer to the current written surgical technique guide, fusionorthopedics.com and associated websites, package insert, product label, and instructions for use before considering and using any Fusion Product. Products and associated accessories may not be available in all markets because of Product availability, and Products are subject to the regulatory law and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Products in your area. Patients must solely rely on advice from their surgeon regarding the risks and all steps of a particular surgery with the utilization of Fusion's Products. As with any surgery, individual results and recovery time vary. Recovery time may range from months to years. Known risks associated with surgery are, but not limited to, injury, significant scarring, nerve damage, mal-union, long-term pain, implant pain, non-union, misalignment, loss of motion, bone shortening, infections, swelling, stiffness, bone fracture, and subsequent required amputation. These works, and its contents, are protected by U.S. and International copyright laws. Reproduction and/or distribution of these works, and its contents, without the express written consent of Fusion Orthopedics USA, LLC is prohibited. © 2023 Fusion Orthopedics USA, LLC