

POLY/LOCK

SMALL BONE PLATES



A TRADEMARK OF
Fusion
ORTHOPEDICS

STRAIGHT PLATES

Thin Profile


Malleable plates provide contoured fit with screws sitting flush on both sides.

Compression Holes

Eccentric holes provide added compression on select plates.

Periosteal Elevator

Leading edge designed to easily remove periosteum from the bone.



COTTON PLATES

Top



Neutral



4mm



6mm




8mm



DUAL MPJ

Universal



WEDGE PLATES

Top



3mm



4mm




5mm



5TH MET PLATES


Jones Plates

Left / Right



Hook Plates

Left / Right



U-PLATES

Left / Right




T-PLATES



Screw Heads

Low profile head prevents stress risers and soft tissue irritation.



Poly Locking

30° poly-axial screw locking for desired angle placement.



EVANS PLATES




Scan to learn about the PolyLock Modular Tray

POLYLOCK

SMALL BONE PLATES

MODULAR TRAY

Instrument Deck




Base Tray

Screw Caddies


SCREWS

2.0mm
8mm - 30mm




Locking
NonLocking

2.4mm
8mm - 30mm



Locking
NonLocking

2.7mm
8mm - 30mm



Locking
NonLocking

4.0mm
18mm - 36mm



CADDIES

Opening Wedge Caddy



Evans Caddy



Small Bone Utility Caddy



5th Met Caddy



Cotton Caddy



FUSION ORTHOPEDICS USA, LLC.
6859 E Rembrandt Ave #122
Mesa, Arizona 85212 USA
(800)-403-6876

Fusion
ORTHOPEDICS

FDA Cleared.
All Rights Reserved.
PolyLock is a Trademark
of Fusion Orthopedics USA, LLC.
SLS016 10032023

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Surgeons should rely on their education and professional training, as the Fusion Instruments, implants, and techniques ("Products") augment traditional surgeries. Fusion Orthopedics, and its affiliates ("Fusion"), does not dispense medical advice and while Fusion recommends that surgeons be trained in the use of any particular Fusion Product before using it in surgery, surgeons must still rely on their professional training and expertise. The information presented is intended to demonstrate the breadth of Fusion's Products. A surgeon and patient must always refer to the current written surgical technique guide, fusionorthopedics.com and associated websites, package insert, product label, and instructions for use before considering and using any Fusion Product. Products and associated accessories may not be available in all markets because of Product availability, and Products are subject to the regulatory law and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Products in your area. Patients must solely rely on advice from their surgeon regarding the risks and all steps of a particular surgery with the utilization of Fusion's Products. As with any surgery, individual results and recovery time vary. Recovery time may range from months to years. Known risks associated with surgery are, but not limited to, injury, significant scarring, nerve damage, mal-union, long-term pain, implant pain, non-union, misalignment, loss of motion, bone shortening, infections, swelling, stiffness, bone fracture, and subsequent required amputation. These works, and its contents, are protected by U.S. and International copyright laws. Reproduction and/or distribution of these works, and its contents, without the express written consent of Fusion Orthopedics USA, LLC is prohibited. © 2023 Fusion Orthopedics USA, LLC