


POLY/LOCK

MIDFOOT PLATES



LISFRANC PLATES

Low Profile
Distal ends of plate are ultra low profile for non-locking screws to reduce soft tissue irritation and to distribute the load.

Thin Profile
Malleable plates provide contoured fit with screws sitting flush

Compression Slot
Self-compressing slots provide added compression on select plates

Large
Medium
Small

MTP PLATES

Raptor Plates
Left / Right

Neutral Dorsal Plates
Left / Right

Universal Plates
Large Medium Small

Chubby Plates
Large Small

Various Hole Sizes
2-8 Holes

2 Hole
3 Hole
4 Hole
5 Hole
6 Hole
8 Hole

LAPIDUS PLATES

Maverick Plates
Left / Right

Medial Lapidus Plates
Left / Right

T-Plates
4 Hole 5 Hole

LAG T-PLATE
Universal

0mm, 2mm - Step Options

L-PLATES

Various Sizes
Plates ranging from 2-11 Hole configurations

4 Hole
5 Hole
6 Hole
6 Hole
5 Hole
4 Hole

Poly Axial Locking
30° Poly-axial screw locking for desired angle placement




Scan to learn about the PolyLock Modular Tray

POLYLOCK

MIDFOOT PLATES

MODULAR TRAY



SCREWS




CADDIES

Utility Caddy



MPJ Caddy



Lisfranc Caddy



Lapidus Caddy



L-Plate Caddy



FUSION ORTHOPEDICS USA, LLC.
6859 E Rembrandt Ave #122
Mesa, Arizona 85212 USA
(800)-403-6876

Fusion
ORTHOPEDICS

FDA Cleared.
All Rights Reserved.
PolyLock is a Trademark
of Fusion Orthopedics USA, LLC.
SLS017 Rev11082023

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Surgeons should rely on their education and professional training, as the Fusion Instruments, implants, and techniques ("Products") augment traditional surgeries. Fusion Orthopedics, and its affiliates ("Fusion"), does not dispense medical advice and while Fusion recommends that surgeons be trained in the use of any particular Fusion Product before using it in surgery, surgeons must still rely on their professional training and expertise. The information presented is intended to demonstrate the breadth of Fusion's Products. A surgeon and patient must always refer to the current written surgical technique guide, fusionorthopedics.com and associated websites, package insert, product label, and instructions for use before considering and using any Fusion Product. Products and associated accessories may not be available in all markets because of Product availability, and Products are subject to the regulatory law and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Products in your area. Patients must solely rely on advice from their surgeon regarding the risks and all steps of a particular surgery with the utilization of Fusion's Products. As with any surgery, individual results and recovery time vary. Recovery time may range from months to years. Known risks associated with surgery are, but not limited to, injury, significant scarring, nerve damage, mal-union, long-term pain, implant pain, non-union, misalignment, loss of motion, bone shortening, infections, swelling, stiffness, bone fracture, and subsequent required amputation. These works, and its contents, are protected by U.S. and International copyright laws. Reproduction and/or distribution of these works, and its contents, without the express written consent of Fusion Orthopedics USA, LLC is prohibited. © 2023 Fusion Orthopedics USA, LLC