

POLY/LOCK

ANKLE PLATES



A TRADEMARK OF
Fusion
ORTHOPEDICS

POSTERIOR FIBULA PLATES


Thickness: Standard = 1.5m

Periosteal Elevators

Leading edge designed to easily remove periosteum from the bone.


Thin Profile

Malleable plates provide contoured fit with screws sitting flush.



Compression Slots

Self-compressing slots provide added compression.




UNIVERSAL LATERAL FIBULA PLATES

Thickness: Standard = 1.5m / *Revision = 2.0m

1/3 TUBULAR UNIVERSAL PLATES

Thickness: Standard = 1.5m / Revision = 2.0m




MEDIAL MALLEOLUS PLATES

Thickness: Standard = 1.5m

POSTERIOR LATERAL TIBIA PLATES

Thickness: Standard = 1.5m




HOOK PLATE

Right / Left



Right / Left




Under the Hood

Low profile (above and below) helps prevent stress risers and soft tissue irritation.



ANATOMICAL LATERAL FIBULA PLATES

Thickness: Standard = 1.5m / Revision = 2.0m



Poly-Locking

30° poly-axial screw locking for desired angle placement.




Scan to learn about the PolyLock Modular Tray

POLYLOCK

ANKLE PLATES

MODULAR TRAY



SCREWS



CADDIES


Posterior Fibula Caddy



Utility Caddy



Posterior/Medial Malleolus Caddy



Lateral Fibula Caddy



FUSION ORTHOPEDICS USA, LLC.
6859 E Rembrandt Ave #122
Mesa, Arizona 85212 USA
(800)-403-6876



FDA Cleared.
All Rights Reserved.
PolyLock is a Trademark
of Fusion Orthopedics USA, LLC.
SLS018 10062023

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Surgeons should rely on their education and professional training, as the Fusion Instruments, implants, and techniques ("Products") augment traditional surgeries. Fusion Orthopedics, and its affiliates ("Fusion"), does not dispense medical advice and while Fusion recommends that surgeons be trained in the use of any particular Fusion Product before using it in surgery, surgeons must still rely on their professional training and expertise. The information presented is intended to demonstrate the breadth of Fusion's Products. A surgeon and patient must always refer to the current written surgical technique guide, fusionorthopedics.com and associated websites, package insert, product label, and instructions for use before considering and using any Fusion Product. Products and associated accessories may not be available in all markets because of Product availability, and Products are subject to the regulatory law and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Products in your area. Patients must solely rely on advice from their surgeon regarding the risks and all steps of a particular surgery with the utilization of Fusion's Products. As with any surgery, individual results and recovery time vary. Recovery time may range from months to years. Known risks associated with surgery are, but not limited to, injury, significant scarring, nerve damage, mal-union, long-term pain, implant pain, non-union, misalignment, loss of motion, bone shortening, infections, swelling, stiffness, bone fracture, and subsequent required amputation. These works, and its contents, are protected by U.S. and International copyright laws. Reproduction and/or distribution of these works, and its contents, without the express written consent of Fusion Orthopedics USA, LLC is prohibited. © 2023 Fusion Orthopedics USA, LLC