

Sgarlato
Silastic System

GAIT



Sgarlato Silastic System

GAIT

The Proven Solution
for Small Joint
Arthroplasty

25

Years of
Successful
Clinical
Outcomes


- **Joint preserving** Keller Arthroplasty technique maintains cartilage on the metatarsal head
- **Facilitates flexion** while maintaining a high strength hinge in the sagittal plane with proven spacer design
- **Ultra High performance**, medical-grade, high-consistency silicone elastomer
- **Sterile packaged** and tray based implants and instruments for ease of use



Hallux Rigidus



Arthroplasty
with GAIT



IMPLANT ORDERING

Product Description	Part Number	Spacer Width	Spacer Length
GAIT Flexible Implant Small	30-100-01	6.0mm	14.5mm
GAIT Flexible Implant Medium	30-100-02	6.6mm	16.0mm
GAIT Flexible Implant Large	30-100-03	7.2mm	17.5mm

FUSION ORTHOPEDICS, LLC.

4135 S. Power Rd. Suite 110
Mesa, Arizona 85212 USA
(800)-403-6876
Sales@fusionorthopedics.com

Fusion
ORTHOPEDICS

FDA Cleared

All Rights Reserved. Sgarlato is
a Trademark of Fusion Orthopedics, LLC.
D-4360-SLS-2 Rev A

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Fusion Orthopedics does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Fusion's product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Fusion Orthopedics product.